

## E. coli Transformation

- 1) Grow cells in LB or 2YT to OD<sub>550</sub> 0.3-0.6.
- 2) Spin 5 minutes at 4,000 RPM and 4° C.
- 3) Resuspend cells on ice in 1/20 volume:  
50 ml culture  
→ 2.5 ml  
store as 500  $\lambda$  aliquots  
LB broth, pH 6.1  
10 % PEG 3350  
5 % DMSO  
10 mM MgCl<sub>2</sub>  
10 mM MgSO<sub>4</sub>
- 4) Ice 10 minutes.
- 5) Mix DNA in 100  $\mu$ l KCM (100 mM KCl, 30 mM CaCl<sub>2</sub>, 50 mM MgCl<sub>2</sub>).  
Keep on ice.
- 6) Add 100  $\mu$ l competent cells to DNA/KCM. (divide in 1/2)  
Mix gently.  
Ice 10-20 minutes.
- 7) Heat shock at room temperature 10 minutes. ~~30 min~~
- 8) Express as Usual. or. 50C at 37° shake 1 hr.  
For amp<sup>r</sup> plasmids I leave at room temp 15 minutes, then plate.  
For other drugs, add 1 ml 2YT and roll at 37 ° C for 1 hour.

keeps indefinitely in the fridge.

To Freeze Competent Cells: Make 10% glycerol at step 4, freeze in dry ice/ethanol bath, and store at -80° C. Thaw cells on ice, then continue with step 5.

Derived from Chung and Miller, 1988, NAR, 16#8, p.3580.

Procedure used at Genentech

I freeze cells without adding glycerol and they work well